



COVID 19 RT PCR DETECTION KIT FOR
SARS-COV-2 **COMMON** STRAIN,
SARS-COV-2 **UK**, **SOUTH AFRICAN** & **P2 BRAZILIAN** VARIANT

4 TESTS – 1 TUBE – 45' PCR

The kit can detect the presence of SARS-COV-2 COVID -19 with RT PCR technology:

Each mix contains specific markers for:

2 Specific genes: Orf 1ab and N gene.

Internal control: genomic.

UK Variant.


SA Variant.

P2 Brazilian Variant.

Number of channels used: 4 (FAM, Cy5, VIC, ROX), for an optimal PCR machine compatibility.

Ability to simultaneously detect the presence of the 4 variants (Common, UK, SA, P2 Br).

"One Step" Workflow



Reagents preparation



Sample preparation



RNA

Load sample



Put on the thermocycler

Contacts

OaCP IE LTD

Phoenix House, Monahan Road,
T12H1XY Cork (IE)
Tel. (IE) +353 (0) 212 376 197
Tel. (IT): +39 051 021 8095

ceo@oacp.it
ask@oacp.it